

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
11 March 2004 (11.03.2004)

PCT

(10) International Publication Number
WO 2004/020673 A1

- | | |
|---|--|
| <p>(51) International Patent Classification⁷:
C12N 15/10</p> <p>(21) International Application Number:
PCT/US2003/026557</p> <p>(22) International Filing Date: 25 August 2003 (25.08.2003)</p> <p>(25) Filing Language:
English</p> <p>(26) Publication Language:
English</p> <p>(30) Priority Data:
60/406,654 28 August 2002 (28.08.2002) US</p> <p>(71) Applicant (for all designated States except US): MILLIPORE CORPORATION [US/US]; 290 Concord Road, Billerica, MA 01821 (US).</p> <p>(72) Inventors; and
(75) Inventors/Applicants (for US only): GABRIELS, Joseph, E., Jr. [US/US]; 41 Cleveland Street, Arlington,</p> | <p>MA 02474 (US). MABUCHI, Masaharu [JP/US]; 29 Lexington Drive, Beverly, MA 01915 (US).</p> <p>(74) Agent: DE LUNA, Renato, M.; Millipore Corporation, 290 Concord Road, Billerica, MA 01821 (US).</p> <p>(81) Designated States (national): AU, CA, JP, KR, US.</p> <p>(84) Designated States (regional): European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR).</p> <p>Declaration under Rule 4.17:
— of inventorship (Rule 4.17(iv)) for US only</p> <p>Published:
— with international search report
— with amended claims</p> <p><i>For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.</i></p> |
|---|--|

(54) Title: COMPOSITIONS OF SOLUTION FOR SEQUENCING REACTION CLEAN-UP

(57) **Abstract:** Wash solution and method for purifying sequencing reaction product. The wash solution comprises an effective amount of guanidine in a low ionic solution to reduce or eliminate the presence of dye terminators in a sequencing reaction product. In its method aspects, the present invention comprises the addition of the wash solution to the sequencing reaction product prior to filtration, followed by filtration to reduce or eliminate unincorporated dye terminators. The purified sequencing products can then be resuspended and transferred to an appropriate substrate for sequencing or further preparation. Dye blobs formed from unincorporated dye terminators no longer interfere with the electropherograms generated upon electrophoresis of the sample.

WO 2004/020673 A1